

HAIYE



HY-BYH series Electronic Nasopharyngo Laryngoscope

HY-BYH1, HY-BYH2, HY-BYH3, HY-BYH4



Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 2 sizes monitors are optional, 3.5 inch and 13.3 inch are touchable screen. HD3b model image processor can be connected to high-definition medical displays for better display effects.; The products are divided into 4 specifications: HY-BYH1, HY-BYH2, HY-BYH3, HY-BYH4



3.5 inch
touch screen
image processor



13.3 inch
touch screen
image processor



HD3b
image processor

Parameter

Field angle	120°, tolerance ±10%
Field depth	3mm~150mm
Rotation	up 0° to 130°, down 0° to 130°, tolerance ±10%
Working time	≥240 mins
Instrument channel	1.2mm~2.6mm (≥Ø1.0mm, HY-BYH1 not compatible)
Length of insert part	360mm, tolerance ±5%



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com