


# HAIYE



## HY-MN2.5 Electronic Ureteroscope



### Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 2 sizes monitors are optional, 3.5 inch and 13.3 inch are touchable screen. The HD3b model image processor can be connected to high-definition medical displays for better display effects.



3.5 inch  
touch screen  
image processor



13.3 inch  
touch screen  
image processor



HD3b  
image processor

### Parameter

Field angle	120°, tolerance ±10%
Field depth	3mm~150mm
Rotation	up / down 275° , tolerance 10% ±
OD of insert part	2.5mm
ID of insert part	1.2mm
Length of insert part	700mm, tolerance ±5%



**Shenzhen Haiye Medical Technology Co.,Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com