

HAIYE



## HY-SCSTX1 Endoscopic Image Processor



### Application

This product consists of an image processing section and is used in conjunction with our company's series of reusable endoscopes for internal surgical area image acquisition, and storage processing during endoscopic surgery.

This product is used in conjunction with the electronic endoscope produced by our company, to process the images captured by the electronic endoscope.

### Feature

This product has complete functions, stable performance, and powerful image processing capabilities; Such as photography, video recording, storage, playback, white balance, etc. It is easy to carry and can seamlessly connect to our reusable endoscope.

### Parameter

Built-in display	3.5", touch screen, 640*480
Compatible models	Reusable bronchoscope, Reusable nasopharyngoscope , Reusable ureteroscope
Video output	HDMI, Type-C
Function	Take pictures, record, freeze, play back, white balance, etc
Language	Optional (English, Spanish, French, Arabic, Russian, Chinese, German, Portuguese, Japanese, Thai, Vietnamese, Korean)
Aspect ratio	4:3
Working temperature	-10°C~40°C
Working power	Voltage: a.c.220V; Frequency: 50Hz; Input power: 30VA
Memory function	Standard with 32GB memory card, storage format:jpg,mov.



**Shenzhen Haiye Medical Technology Co.,Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com