



HYHJ-1330 Disposable Video laryngoscope

/HYHJ-1330 I, HYHJ-1330 II, HYHJ-1330 III/



Application

It is provided for clinical use to provoke the patient's epiglottis to expose the glottis, and to guide medical staff to perform tracheal intubation and implement mechanical ventilation.

Suitable for anesthesia departments, ICU, pediatrics, emergency centers, ambulances, etc. in medical institutions.



Feature and Configuration

This product has a simple and durable structure, and is designed as an all-in-one machine with one specification, applicable to Adults, children, newborns; corresponding to 3 kinds of disposable laryngoscope blade, can be applied to all people.



Parameter

Monitor	3" color display
Camera	Pixel: ≥ 2 million; illuminance: ≥ 3000 lux
Battery	Battery life: ≥ 6 h; charging time: ≤ 4 h; capacity: 3400mAh
Function	Take pictures, video, freeze, store and export data
Storage	-40°C~70°C ambient temperature
Humidity	$\leq 93\%$ rh



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com