

HAIYE



HYHJ-KM **Disposable** **Video Laryngoscope** /HYHJ-KM3C, HYHJ-KM5C/



Application

It is provided for clinical use to provoke the patient's epiglottis to expose the glottis, and to guide medical staff to perform tracheal intubation and implement mechanical ventilation.

Suitable for anesthesia departments, ICU, pediatrics, emergency centers, ambulances,



Feature and Configuration

This product is divided into two models: HYHJ-KM3C and HYHJ-KM5C, has a total of 8 types of disposable laryngoscope blades (Miller00, Miller0, Miller1, Mac1, Mac2, Mac3, Mac4, Mac5), can be applied to all people.

Single-use laryngoscope blade



Parameter

Monitor	3" LCD display / 3.5" color touch screen display
Camera	pixel:≥3million; light illuminance:≥3000 lux
Battery	battery life:≥4h;charging time:≤5h;capacity:3400mAh
Interface Function	HDMI,Type-C;Photo,video,freeze,playback,store and export data
Storage	-40°C~70°C ambient temperature
Humidity	≤93% rh



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com