

HAIYE

HY-ZQ series Electronic Bronchoscope

/HY-ZQ2.0, HY-ZQ2.8, HY-ZQ4.5, HY-ZQ5.2, HY-ZQ5.8/



Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 2 sizes monitors are optional, 3.5 inch and 13.3 inch touch screen. The HD3b model image processor can be connected to high-definition medical displays for better display effects.

The products are divided into 5 specifications: HY-ZQ2.0, HY-ZQ2.8, HY-ZQ4.5, HY-ZQ5.2 and HY-ZQ5.8.



3.5 inch
touch screen
image processor



13.3 inch
touch screen
image processor



HD3b
image processor

Parameter

Field angle	120°, tolerance $\pm 10\%$
Field depth	3mm~150mm
Rotation	up 0° to 180°, down 0° to 130°, tolerance $\pm 10\%$
Working time	≥ 240 mins
Instrument channel	1.2mm~3.0mm ($\geq \varnothing 1.0$ mm, HY-ZQ2.0 not compatible)
Length of insert part	600mm, tolerance $\pm 5\%$



Shenzhen Haiye Medical Technology Co., Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com