

# HAIYE



## HY-TC1 Mobile Stand




### Application

Equipped with electronic endoscope image processor produced by our company, it can be used with our video laryngoscope, disposable uroscope, bronchoscope and other products for clinical tracheal intubation or treatment.

It is suitable for ENT, operating room, ICU, pediatrics, emergency center, ambulance, respiratory medicine, infection department, thoracic surgery, urology and so on.

### Available Product



### Product Description

The Mobile Stand consists of a support plate, support rods, workbench, multi-layer racks, bottom storage box, chassis, and universal casters. It can carry display screens of various specifications. The rack can hold image processors, endoscope camera systems, endoscope cold light sources, and printers, etc.



**Shenzhen Haiye Medical Technology Co., Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: [info@haiyemedical.com](mailto:info@haiyemedical.com)

Website: [www.hyendoscope.com](http://www.hyendoscope.com)