

HAIYE



HY-YM2.8

Single-use Flexible Electronic Ureteroscope

HY-YM2.8A, HY-YM2.8B



Application

Used in conjunction with the electronic endoscope image processor produced by our company, it can be inserted into the urethra and renal pelvis through the urethra for observation, diagnosis, photography and treatment of the urethra and renal pelvis.



Configuration

The product is a disposable handheld flexible electronic endoscope consisting of an electronic ureteropyeloscope, a perfusion plug and an optional accessory waterproof cover.



Video converter



13.3 inch touch screen image processor



HD3b image processor

Parameter

Field angle	120°, tolerance ±15%
Field depth	3 ~100mm
Rotation	up/down 275°, tolerance ±15%
OD of insert part	2.5mm
ID of insert part	≥1.2mm
Length of insert part	700mm



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com