

# HAIYE



## HY-YPG

### Single-use Flexible Electronic Cystoscope

HY-YPG3.9, HY-YPG-4.5, HY-YPG5.3



### Application

Used in conjunction with the electronic endoscope image processor produced by our company. It is inserted into the urethra and renal pelvis through the urethra and is used for observation, diagnosis, photography and treatment of the urethra and renal pelvis.



### Configuration

This product is a disposable handheld flexible electronic endoscope. The operating part is designed in accordance with ergonomics. It can take pictures, record videos, and switch freely. The whole machine is a disposable design, which effectively prevents cross infection.



Video converter



13.3 inch touch screen image processor



HD3b image processor

Field angle	120°, tolerance ±10%
Field depth	3~100 mm
Camera pixels	400*400
Product specifications	HY-YPG3.9; HY-YPG4.5; HY-YPG5.3
ID of insert part	≥1.4-2.8 mm ± 0.2 mm
OD of insert part	≥3.9-5.3 mm ± 0.3 mm



**Shenzhen Haiye Medical Technology Co.,Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: [info@haiyemedical.com](mailto:info@haiyemedical.com)

Website: [www.hyendoscope.com](http://www.hyendoscope.com)