

HAIYE



HYHJ-KX Series

Disposable Visual Laryngoscope Blade



Models:

HYHJ-KX00a, HYHJ-KX00b
HYHJ-KX0a, HYHJ-KX0b
HYHJ-KX1a, HYHJ-KX1b
HYHJ-KX2a, HYHJ-KX2b
HYHJ-KX3a, HYHJ-KX3b
HYHJ-KX4a, HYHJ-KX4b

Application

For clinical provocation of the patient's epiglottis to expose the vocal folds, guiding healthcare professionals to extubate and implement mechanical ventilation under visualization.

Applicable to the anesthesiology department, ICU, pediatrics department of medical institutions, emergency center, ambulance, etc.



Feature

The product is aseptically packaged for single use, minimizing the risk of cross infection.

The product is available in 6 sizes, and work with our video processors, it can be used by most people, from pediatric to adult.



Separate cable
and laryngoscope blade



one piece
(cable and laryngoscope blade)

Parameter

Monitor	12.1" or 13.3" touch screen (purchase separately)
Camera	one million pixels
Function	Take picture, record, freeze, playback, store and export data
Storage	ambient environment -40°C~70°C
Humidity	≤93% rh



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com