

HAIYE



HY-YBYH series Single-use Nasopharyngoscope

/HY-YBYH2.0, HY-YBYH2.5, HY-YBYH3.9, HY-YBYH4.5, HY-YBYH5.3 /



Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 2 sizes monitors are optional: video converter and 13.3 inch are touchable screen. The HD3b model image processor can be connected to high-definition medical displays for better display effects.

The products are divided into 5 specifications:

HY-YBYH2.0, HY-YBYH2.5, HY-YBYH3.9, HY-YBYH4.5, HY-YBYH5.3.



Video converter



13.3 inch
touch screen
image processor



HD3b
image processor

Parameter

| | |
|-----------------------|--|
| Field angle | 120°, tolerance ±10% |
| Field depth | 3mm~100mm |
| Rotation | up 0° to 210°, down 0° to 210, tolerance ±10% |
| Working time | ≥240 mins |
| Instrument channel | 1.2mm~2.8mm (≥Ø1.0mm, HY-YBYH2.0 not compatible) |
| Length of insert part | 380mm ,tolerance ±5% |



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com