

HAIYE



HY-YZQ series

Disposable Bronchoscope

/HY-YZQ2.0, HY-YZQ2.5, HY-YZQ4.5, HY-YZQ5.3/



Usage

Used in conjunction with the electronic endoscope image processor produced by our company, it is used by medical units to provide the necessary images for observation, diagnosis, photography and bronchi and treatment of the trachea and bronchi of adults and children over 3 months old. Product structure: The product is mainly composed of bending part, inserting part, operation part, conical sheath and video connection line.



Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 2 sizes monitors are optional: video converter and 13.3 inch touch screen. The HD3b model image processor can be connected to high-definition medical displays for better display effects.

The products are divided into 4 specifications: HY-YZQ2.0, HY-YZQ2.5, HY-YZQ4.5 and HY-YZQ5.3.



Video converter



13.3 inch touch screen image processor



HD3b image processor

Parameter

Field angle	120°, tolerance ±10%
Field depth	3 ~100mm
Rotation	up 0° to 180°, down 0° to 130°
Camera	CMOS electronic imaging
Instrument channel	1.2mm~ 2.8mm (≥Ø1.0mm, HY-YZQ2.0 not compatible)
Length of insert part	600mm, tolerance ±5%



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com