

# HAIYE



## HY-RD series Electronic Choledochoscope

HY-RD2.0 , HY-RD2.8 , HY-RD4.5 , HY- RD5.2



### Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 2 sizes monitors are optional, 3.5 inch and 13.3 inch are touchable screen. HD3b model image processor can be connected to high-definition medical displays for better display effects.; The products are divided into 4 specifications: HY-RD2.0 , HY-RD2.8 , HY-RD4.5 , HY- RD5.2



3.5 inch  
touch screen  
image processor



13.3 inch  
touch screen  
image processor



HD3b  
image processor

### Parameter

Field angle	120°, tolerance ±10%
Feild depth	3mm~150mm
Rotation	up 0° to 130°, down 0° to 130°, tolerance ±10%
Working time	≥240 mins
Instrument channel	1.2mm~2.6mm (≥Ø1.0mm, HY-BYH1 not ompatible)
Length of insert part	360mm, tolerance ±5%



**Shenzhen Haiye Medical Technology Co.,Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com